



# Aspiration Containment System for Large Ovarian Cyst Removal

Atiya Anwar, Christabel Ekeocha, Emily Graba, Mary William, Kayla Whatley

Clinical Advisor: Dr. Jordan Mattson, Industry Advisor: Joram Slager




## Need Statement

Surgeons require a simple solution for aspirating large ovarian cysts via a small surgical incision and without any content leakage. Such a solution will significantly reduce post-operative recovery times and risk of cancer metastasis.

## Background

### Clinical Problem:

- Nearly 21,000 women in the United States develop a cancerous ovarian cyst annually, resulting in 14,500 deaths a year.
- Ovarian cysts can grow to be upwards of 30L in volume.
- There is no reliable method to differentiate between a malignant and benign ovarian cyst in vivo.
- It is critical to remove cysts without leakage of content into the abdomen, thereby eliminating the potential for metastasis and need for chemotherapy.



<https://www.starambhartia.org/blog/gynecology/ovarian-cystectomy/>

**Current Technology:** There are two main surgical techniques:


1. **Laparotomy:** An incision is made from the pubis to the sternum and the cyst is removed whole from the body (i.e. without puncture).
2. **Laparoscopy:** A small incision is made and a containment bag is maneuvered around the cyst. The cyst is aspirated in the open bag using the Neptune 3 Waste Management System.

Extensive post-op recovery times

Technical challenges and high risk of leakage

## Prototyping

### Timeline:




**Process:** Several iterations were created between prototypes:

- Single and dual chamber systems
- With mesh (bottom) and containment bag (top)
- Varying device heights
- Varying locations of vacuum connection post

Prototypes were initially made from scratch materials, then later 3D printed with PLA, ABS, and finally PolyJet material.

## Device Design



The ACS utilizes negative pressure to create an airtight seal with the surface of a cyst. The device is composed of 2 chambers:

- **Inner Chamber:** allows the surgeon to place and utilize aspiration tools they are well acquainted with (e.g. scalpel & yankauer suction tip) in a comfortable-to-manuever space.
- **Outer Chamber:** hosts a pathway for negative pressure to be delivered from a surgical vacuum to the device base.

The soft bottom ensures a smooth interface that promotes cyst adhesion. The septum serves as an area for which sterile buffer solution can be injected to clean the outer chamber.

## Preliminary Results


### Cyst Modeling:

- Initially, water-filled balloons and plastic Ziploc bags were used.
- Ballistic Gel #1 from Humimic Medical was molded into a hollow spherical sack, then filled with red-dyed water.
- VOC feedback indicated that ballistic gel realistically mimics cyst material properties!



### Device Testing:

- **Leak Tests**
  - No leakage of fluid occurred during or after aspiration.
  - Device creates a strong, airtight seal with cyst surface.
- **Compressive Force Tests**
  - Device body can withstand 80 lbf (avg. hand grip = 60 lbf).
  - Vacuum post can withstand 35 lbf.
- **VOC Tests**
  - Device is intuitive to use, compatible with multiple removal techniques, and provides adequate visualization.



## Future Research Directions

1. Incorporation of a pressure release valve
2. Selection of an autoclavable, non-cytotoxic, and biocompatible device material
3. Modification of the vacuum tubing for compatibility with the Neptune 3 Waste Management System
4. Further evaluation of device ergonomics

## Acknowledgments

The team would like to thank the following people:

- Cyrus Rezvanifar
- Dr. Jordan Mattson
- Joram Slager
- Dr. Peter Argenta
- Dr. Andrea O'Shea
- Ben Guengerich
- Anderson Labs staff
- Protolabs
- Mark Sanders
- UMN MVP Challenge